

## Hot+Cold therapy Kinetec Kooler™ Consumables

The Kinetec Kooler™ range helps to relieve pain and reduce swelling through a wide variety of carefully selected items providing hot and cold therapy. Use for minor injuries, joint pain, muscle aches and pain.



### Kinetec Kooler™ hot/cold gel pack

This thermotherapy pack can be used to heat or cool all parts of the body. Featuring two adjustable fixing straps that can be used to secure the pack in place. Pack size is 13,5x23,5cm. The pack is filled with non-toxic gel, does not drip and remains flexible. The pack is reusable and is easy to clean and decontaminate. The temperature of the pack can be controlled by immersion in water or by using a microwave or a freezer. The individual box size is 10x5,5x14cm.

Ordering information: 4665009807 (box of 1)

### Kinetec Kooler™ instant cold pack

The instant cold pack is single use. Activation is immediate and performed by squeezing the pack as indicated. Pack size is 12x23cm. Size of the 24-item box is 23x23x22cm.

Ordering information: 4665009815 (box of 24)



### Kinetec Kooler™ instant hot pack

The Instant hot pack is single use, it provides instantaneous heat up to 43°C (110°F). It remains effective for approximately 20 to 25 minutes. Usable anywhere and at any time. Activation is immediate and performed by squeezing the pack. Contents – non-toxic: Sodium Thiosulfate, water.

Pack size is 12x23cm.

Size of 24-item box is 22x13,5x22cm.

Ordering information: 4665010028 (box of 24)



### Kinetec Kooler™ Arktic Breeze spray

Quickly cools down the skin and keeps it cool to relieve pain and reduce swelling. Works in all directions. Size of the spray is 18x4,5cm, 150ml.

Ordering information: 4665009823 languages DE and FR (box of 12)

4665009774 languages EN and NL (box of 12)

4665009831 languages IT and ES (box of 12)



### Kinetec Kooler™ Arktic Breeze gel sachet

Quickly cools down the skin and keeps it cool to relieve pain and reduce swelling.

Size of the sachet is 8x10cm, 5ml. Size of the 30-item box is 16x10x5,5cm.

Ordering information: 4665009782 languages EN - NL and ES (box of 30)

4665009790 languages DE - FR and IT (box of 30)



### Kinetec Kooler™ Arktic Breeze gel tube

Quickly cools down the skin and keeps it cool to relieve pain and reduce swelling. Size of the individual box is 5,5x5x15,5cm, 100ml.

Ordering information: 4665009849 languages EN and NL (box of 12)

4665009857 languages DE and FR (box of 12)

4665009865 languages ES and IT (box of 12)

